# Home Hemodialysis with the Tablo System: The First 1,000 Real-World Treatments

Michael Aragon, MD Josh Schumacher

# **.**Outset

# **BACKGROUND:**

The Tablo<sup>®</sup> Hemodialysis System is an all in one, easy-to-learn device capable of achieving clearance goals in as little as three treatments per week. It features integrated water purification, on-demand dialysate production and two-way wireless data transmission.

Tablo obtained FDA clearance for Home Hemodialysis (HHD) in March 2020 based on a prospective, cross-over between in-center and at home dialysis trial (NCT02460263).

Tablo met all safety and effectiveness endpoints, reported high rates of treatment adherence and patient retention, with greatly reduced training time.

# **OBJECTIVE:**

Report on the first 1,000 treatments performed on Tablo by patients at home in the real-world.

#### **METHODS:**

- Incident and prevalent patients currently performing in-center dialysis, PD, or HHD were initiated on Tablo.
- Patients underwent training by facility nurses prior to beginning treatments at home.
- Data on the first 1000 treatments was obtained wirelessly via Tablo's data platform, along with corresponding patient training data.

# **RESULTS:**

- Average training time was 7.4 training days.
- Patient retention was 100% with no patients opting out of HHD with Tablo.
- Mean prescribed treatment time was 3.2 hours with a mean frequency of 3.7 treatments per week.
- Mean UF rate per treatment was 7.3 mL/kg/hr.
- Treatment adherence was 93%, with 95% of treatments completing within 10% of prescribed time.
- The mean number of clinically significant alarms per treatment was 1.0 (± 3.0), with an average time to resolution of 10.7 (± 18.5) seconds.

# **Treatment Results**

| Variable                    | Treatment Data<br>(Mean ± SD) |
|-----------------------------|-------------------------------|
| Patient Weight (kg)         | 89.8 ± 25.4                   |
| Actual Treatment Time (hrs) | 3.2 ± 0.8                     |
| Actual UF Rate (mL/kg/hr)   | 7.3 ± 4.1                     |

This data supports that Tablo is capable of successfully achieving clinical goals while reducing the overall patient burden often associated with HHD.

# **CONCLUSION:**

Results from the Tablo IDE trial are reproducible in the real world at a frequency of 3-4 treatments/week.



- 27.4 hours over 7.9 sessions of Tablo Training
- ✓ 95.4% Treatment Success

✓ 100% Retention Rate

**Prescribed HHD at** 

3.7 treatments/week

- ✓ 93.2% Treatment Adherence
- ✓ 1.0 Alarm per Treatment
- ✓ 10.7 seconds to Alarm Resolution

Observed averages for first 20 patients and 1,000 treatments.